

Просеивающие машины орбитальные, лабораторные сита, тестовые сита, крышки и поддоны для сит

Технические характеристики

По вопросам продаж и поддержки обращайтесь:

Алматы (727)345-47-04
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Саранск (8342)22-96-24
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(727)345-47-04

Беларусь +(375)257-127-884

Узбекистан +998(71)205-18-59

Киргизия +996(312)96-26-47

эл.почта: cen@nt-rt.ru || сайт: <http://coleparmer.nt-rt.ru/>

Cole-Parmer® Orbital Sieve Shakers



Back-and-forth lateral motion combined with up/down and tilting actions distribute material across the sieve

- Accurately separate particles from size no. 4 to no. 200 (4.75 mm to 75 μ m)
- Precisely time operation up to 99 minutes with the built-in digital countdown timer
- EZ-Clamp system quickly and securely locks in sieves for testing
- Compliant with multiple ASTM and AASHTO standards

Use the orbital sieve shakers for particle size testing of mineral aggregates and for sieve gradation of free-flowing granular materials. The back-and-forth and up/down motions plus tilting action move the testing samples across the entire mesh area of the test sieve to quickly separate the materials. Orbital shakers accommodate from no. 4 to no. 200 material size range but may be used with larger sizes when specifications permit.

Shakers feature a digital timer with a large LED display that can be seen from a distance. Use the countdown timer for operations up to 99 minutes with 1 second increment. Timer includes a pause function and user-selectable modes for great versatility. Easy one-button operation and saved settings offer excellent repeatability between cycles.

Quickly and easily lock in the sieves using the EZ-Clamp system. The sieve cover slides freely up and down the clamp rods to hold sieves in place. A quick twist of the clamping screws secures the stack in place. All mechanical and electrical parts are enclosed in a rugged painted steel case for safety and durability.

item	Sieves Size Accepted (OD)
 EW-59998-62	8" or 200 mm

item EW-59998-63	Sieves Size Accepted (OD) 8", 12", 200 mm, or 300 mm
---	--

Item # EW-59998-62

Specifications & Description

- Sieves Size Accepted (OD) 8" or 200 mm
- Timer Range (min) Up to 99 minutes
- Width (InchShort) 17
- Depth (in) 11
- Height (in) 30
- Width (cm) 43.2
- Depth (cm) 27.9
- Height (cm) 76.2
- Power (VAC) 115
- Power (Hz) 60
- Shipping Weight (kg) 27.2
- Shipping Weight (lbs) 60
- DescriptionOrbital Sieve Shaker with Digital Timer, Accepts 8" or 200 mm Diameter Sieves; 115 VAC
- Warranty1 Year

Item # EW-59998-63

Specifications & Description

- Sieves Size Accepted (OD) 8", 12", 200 mm, or 300 mm
- Timer Range (min) Up to 99 minutes
- Width (InchShort) 26
- Depth (in) 20
- Height (in) 40
- Width (cm) 66
- Depth (cm) 50.8
- Height (cm) 101.6
- Power (VAC) 115
- Power (Hz) 60
- Shipping Weight (kg) 45.3
- Shipping Weight (lbs) 100
- DescriptionOrbital Sieve Shaker with Digital Timer, Accepts 8", 12", 200 mm, or 300 mm Diameter Sieves; 115 VAC
- Warranty1 Year

Cole-Parmer Sieve Shaker, Accepts up to 7 Full-Height 3" Sieves; 115V



Cole-Parmer – Item # EW-59986-01

Quiet, electromagnetic vibratory action for small samples

- 0 to 100% amplitude control for fine particle separation
- 3 year warranty

Specifications & Description

- Sieves Size Accepted (OD) 3"
- Timer Range (min) 99
- Vibrational Frequency (VPM) 3600
- Frequency 60 taps/min
- Width (InchShort) 8
- Depth (in) 11
- Height (in) 22

- Width (cm) 20.3
- Depth (cm) 27.9
- Height (cm) 55.9
- Power (VAC) 115
- Power (Hz) 50/60
- Shipping Weight (kg) 7.3
- Shipping Weight (lbs) 16
- Description Sieve Shaker, Accepts up to 7 Full-Height 3" Sieves; 115V
- Warranty 3 Years

MORE ABOUT THIS ITEM

This sieve shaker accommodates up to seven full-height 3"-diameter sieves and is ideal for small samples of chemicals, minerals, pharmaceuticals, powdered metals, cosmetics, abrasives, ores, foods and other fine powders. The stainless steel, 3"-diameter sieves, sieve cover, and sieve pan are sold separately.

High frequency, electromagnetic vibratory action with 0 to 100% amplitude control is ideal for fine particle separations. Inline solenoid-actuated tapping action may be used during vibration to speed dry separations. Also use the tap separately for tap-settling and bulk density tests.

The unique clamping mechanism makes insertion and removal of the sieves simple and convenient and holds the sieves securely during operation. Shaker also features on/off controls for vibration and a manual timed-mode switch. Digital timer shows time set and time remaining to enable reliable, repeatable results

Cole-Parmer One-Touch Vibratory Sieve Shaker 115, 60 (VAC, Hz)



Cole-Parmer – Item # EW-59986-03

of respondents would recommend this to a friend

Touch-screen programmable controller allows selection of vibration level, time, and pauses

- Store up to 99 testing profiles for easy repeat of standard tests
- 3 year warranty

Specifications & Description

- Sieves Size Accepted (OD) 8"
- Timer Range (min) 99
- Width (InchShort) 16 1/2
- Depth (in) 12
- Height (in) 28
- Width (cm) 42
- Depth (cm) 31
- Height (cm) 71.1
- Power (VAC) 115
- Power (Hz) 60
- Shipping Weight (kg) 43
- Shipping Weight (lbs) 94
- DescriptionOne-Touch Vibratory Sieve Shaker
- Warranty3 Years



MORE ABOUT THIS ITEM

Combine proven separation technology with the latest in electronic control for fast, accurate separations. Establish the particle size for aggregates, ores, pharmaceuticals, powders, and more. Simply select vibration level, time, and pauses, and press "Start." Three-dimensional sieving action evenly distributes the sample across the sieving surface and reorients particles to ensure the maximum number of passing opportunities. Store up to 99 testing profiles for your standard sieving operations to ensure repeatability.

Shaker accepts eight full-height or sixteen half-height 8" dia. sieves. Sieve clamps allow fast, easy insertion and removal of sieve stacks, while ensuring that sieves are tightly secured during testing.

Cole-Parmer® Tapping Sieve Shakers



Quickly produces sharp, repeatable separations

- Vertical two-way tapping at 48 taps/minute
- Innovative EZ-Clamp system secures sieve stacks
- Internal counterbalance system
- Precision digital controller
- Accommodates up to 20 sieves

These tapping sieve shakers allow you to produce sharp, repeatable separations quickly over a broad range of particle sizes and material types. These high-capacity shakers can hold as many as 20 sieves. The sieving action results from circular rotation of the sieve stack along with vertical two-way tapping at 48 taps/minute,

item	Power (VAC)	Shipping Weight (kg)	Sieves Size Accepted (OD)
 EW-59998-60	115	61.24	8" (200 mm)

item	Power (VAC)	Shipping Weight (kg)	Sieves Size Accepted (OD)
 EW-59998-61	115	73.48	8" & 12" (200 & 300 mm)

Item # EW-59998-61

Specifications & Description

- Sieves Size Accepted (OD) 8" & 12" (200 & 300 mm)
- Timer Range (min) Up to 99 minutes, 59 seconds
- Width (InchShort) 22 1/2
- Depth (in) 13
- Height (in) 41
- Width (cm) 57.2
- Depth (cm) 33
- Height (cm) 104.1
- Power (VAC) 115
- Power (Hz) 60
- Shipping Weight (kg) 73.48
- Shipping Weight (lbs) 162
- DescriptionTapping Sieve Shaker, 8"/12"; 115V
- Warranty1 year

Item # EW-59998-60

Specifications & Description

- Sieves Size Accepted (OD) 8" (200 mm)
- Timer Range (min) Up to 99 minutes, 59 seconds
- Width (InchShort) 22 1/2
- Depth (in) 13
- Height (in) 41
- Width (cm) 57.2
- Depth (cm) 33
- Height (cm) 104.1
- Power (VAC) 115
- Power (Hz) 60
- Shipping Weight (kg) 61.24
- Shipping Weight (lbs) 135
- DescriptionTapping Sieve Shaker, 8"; 115V
- Warranty1 year

Cole-Parmer® Testing Sieves, 12" dia, Stainless Steel Frame/Stainless Steel Wire



Get more accurate results

- Sievers are serialized
 - Come with a test sieve certificate of compliance
- These stainless steel (SS) testing sieves are serialized and come with a test sieve certificate of compliance. Equivalent to ASTM E11-09 compliance sieve.

•

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
EW-59994-00	3.25	2.36	8

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59994-01	3.25	2	10
 EW-59994-02	3.25	1.68	12
 EW-59994-03	3.25	1.41	14

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59994-04	3.25	1.19	16
 EW-59994-05	3.25	1	18
 EW-59994-06	3.25	0.841	20

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59994-07	3.25	0.707	25
 EW-59994-08	3.25	0.595	30
 EW-59994-09	3.25	0.5	35

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59994-10	3.25	0.425	40
 EW-59994-11	3.25	0.355	45
 EW-59994-12	3.25	0.297	50

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59994-13	3.25	0.25	60
 EW-59994-14	3.25	0.212	70
 EW-59994-15	3.25	0.177	80

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59994-16	3.25	0.149	100
 EW-59994-17	3.25	0.125	120
 EW-59994-18	3.25	0.105	140

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59994-19	3.25	0.088	170
 EW-59994-20	3.25	0.074	200
 EW-59994-21	3.25	0.063	230

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59994-22	3.25	0.053	270
 EW-59994-23	3.25	0.044	325
 EW-59994-24	3.25	0.037	400

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59995-00	1.625	2.36	8
 EW-59995-01	1.625	2	10
 EW-59995-02	1.625	1.68	12

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59995-03	1.625	1.41	14
 EW-59995-04	1.625	1.19	16
 EW-59995-05	1.625	1	18

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59995-06	1.625	0.841	20
 EW-59995-07	1.625	0.707	25
 EW-59995-08	1.625	0.595	30

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59995-09	1.625	0.5	35
 EW-59995-10	1.625	0.425	40
 EW-59995-11	1.625	0.355	45

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59995-12	1.625	0.297	50
 EW-59995-13	1.625	0.25	60
 EW-59995-14	1.625	0.212	70

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59995-15	1.625	0.177	80
 EW-59995-16	1.625	0.149	100
 EW-59995-17	1.625	0.125	120

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59995-18	1.625	0.105	140
 EW-59995-19	1.625	0.088	170
 EW-59995-20	1.625	0.074	200

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59995-21	1.625	0.063	230
 EW-59995-22	1.625	0.053	270
 EW-59995-23	1.625	0.044	325

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59995-24	1.625	0.037	400

Cole-Parmer® Testing Sieves, 8" dia, Stainless Steel Frame/Stainless Steel Wire



Don't second guess accuracy

- Come with test sieve certificate of compliance
 - Serialized
- These stainless steel (SS) testing sieves are serialized and come with a test sieve certificate of compliance. Equivalent to ASTM E11-09 compliance sieve.

•

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
EW-59984-00	2	2.36	8

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59984-01	2	2	10
 EW-59984-02	2	1.68	12
 EW-59984-03	2	1.41	14

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
A cylindrical metal sieve with a fine mesh, shown from a top-down perspective. EW-59984-04	2	1.19	16
A cylindrical metal sieve with a fine mesh, shown from a top-down perspective. EW-59984-05	2	1	18
A cylindrical metal sieve with a fine mesh, shown from a top-down perspective. EW-59984-06	2	0.841	20

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59984-07	2	0.707	25
 EW-59984-08	2	0.595	30
 EW-59984-09	2	0.5	35

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59984-10	2	0.425	40
 EW-59984-11	2	0.355	45
 EW-59984-12	2	0.297	50

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59984-13	2	0.25	60
 EW-59984-14	2	0.212	70
 EW-59984-15	2	0.177	80

item		Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59984-16		2	0.149	100
 EW-59984-17		2	0.125	120
 EW-59984-18		2	0.106	140

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59984-19	2	0.088	170
 EW-59984-20	2	0.074	200
 EW-59984-21	2	0.063	230

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59984-22	2	0.053	270
 EW-59984-23	2	0.044	325
 EW-59984-24	2	0.037	400

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
A cylindrical metal sieve with a fine mesh. EW-59985-00	1	2.36	8
A cylindrical metal sieve with a fine mesh. EW-59985-01	1	2	10
A cylindrical metal sieve with a fine mesh. EW-59985-02	1	1.68	12

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59985-03	1	1.41	14
 EW-59985-04	1	1.19	16
 EW-59985-05	1	1	18

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59985-06	1	0.841	20
 EW-59985-07	1	0.707	25
 EW-59985-08	1	0.595	30

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59985-09	1	0.5	35
 EW-59985-10	1	0.425	40
 EW-59985-11	1	0.355	45

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59985-12	1	0.297	50
 EW-59985-13	1	0.25	60
 EW-59985-14	1	0.212	70

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59985-15	1	0.177	80
 EW-59985-16	1	0.149	100
 EW-59985-17	1	0.125	120

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59985-18	1	0.105	140
 EW-59985-19	1	0.088	170
 EW-59985-20	1	0.074	200

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59985-21	1	0.063	230
 EW-59985-22	1	0.053	270
 EW-59985-23	1	0.044	325

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59985-24	1	0.037	400

Cole-Parmer® Testing Sieves, 8" dia, Brass Frame/Stainless Steel Wire



- Economical brass frame combined with the durability of stainless steel mesh
These brass and stainless steel (SS) testing sieves are serialized and come with a test sieve certificate of compliance. Equivalent to ASTM E11-09 compliance sieve.

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
EW-59940-00	2	2.36	8

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59940-01	2	2	10
 EW-59940-02	2	1.68	12
 EW-59940-03	2	1.41	14

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59940-04	2	1.19	16
 EW-59940-05	2	1	18
 EW-59940-06	2	0.841	20

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59940-07	2	0.707	25
 EW-59940-08	2	0.595	30
 EW-59940-09	2	0.5	35

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59940-10	2	0.425	40
 EW-59940-11	2	0.355	45
 EW-59940-12	2	0.297	50

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59940-13	2	0.25	60
 EW-59940-14	2	0.212	70
 EW-59940-15	2	0.177	80

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59940-16	2	0.149	100
 EW-59940-17	2	0.125	120
 EW-59940-18	2	0.105	140

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59940-19	2	0.088	170
 EW-59940-20	2	0.074	200
 EW-59940-21	2	0.063	230

item		Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59940-22		2	0.053	270
 EW-59940-23		2	0.044	325
 EW-59940-24		2	0.037	400

item		Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59941-00		1	2.36	8
 EW-59941-01		1	2	10
 EW-59941-02		1	1.68	12

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59941-03	1	1.41	14
 EW-59941-04	1	1.19	16
 EW-59941-05	1	1	18

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59941-06	1	0.841	20
 EW-59941-07	1	0.707	25
 EW-59941-08	1	0.595	30

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59941-10	1	0.425	40
 EW-59941-12	1	0.297	50
 EW-59941-13	1	0.25	60

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59941-14	1	0.212	70
 EW-59941-15	1	0.177	80
 EW-59941-16	1	0.149	100

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59941-17	1	0.125	120
 EW-59941-18	1	0.105	140
 EW-59941-20	1	0.074	200

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59941-22	1	0.053	270
 EW-59941-23	1	0.044	325
 EW-59941-24	1	0.037	400

Cole-Parmer® Testing Sieves, 8" dia, Brass Frame/Brass Wire



Get more accurate results

- Sievers are serialized
 - Come with a test sieve certificate of compliance
- These brass testing sieves are serialized and come with a test sieve certificate of compliance. Equivalent to ASTM E11-09 compliance sieve.
-

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
EW-59980-00	2	2.36	8

item		Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59980-01		2	2	10
 EW-59980-02		2	1.68	12
 EW-59980-03		2	1.41	14

item		Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59980-04		2	1.19	16
 EW-59980-05		2	1	18
 EW-59980-06		2	0.841	20

item		Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59980-07		2	0.707	25
 EW-59980-08		2	0.595	30
 EW-59980-09		2	0.5	35

item		Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59980-10		2	0.425	40
 EW-59980-11		2	0.355	45
 EW-59980-12		2	0.297	50

item		Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59980-13		2	0.25	60
 EW-59980-14		2	0.212	70
 EW-59980-15		2	0.177	80

item		Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59980-16		2	0.149	100
 EW-59980-17		2	0.125	120
 EW-59980-18		2	0.105	140

item		Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59980-19		2	0.088	170
 EW-59980-20		2	0.074	200
 EW-59980-21		2	0.063	230

item		Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59980-22		2	0.053	270
 EW-59980-23		2	0.044	325
 EW-59981-00		1	2.36	8

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59981-01	1	2	10
 EW-59981-02	1	1.68	12
 EW-59981-03	1	1.41	14

item		Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59981-04		1	1.19	16
 EW-59981-05		1	1	18
 EW-59981-06		1	0.841	20

item		Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59981-07		1	0.707	25
 EW-59981-08		1	0.595	30
 EW-59981-09		1	0.5	35

item		Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59981-10		1	0.425	40
 EW-59981-11		1	0.355	45
 EW-59981-12		1	0.297	50

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59981-13	1	0.25	60
 EW-59981-14	1	0.212	70
 EW-59981-15	1	0.177	80

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59981-16	1	0.149	100
 EW-59981-17	1	0.125	120
 EW-59981-18	1	0.105	140

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59981-19	1	0.088	170
 EW-59981-20	1	0.074	200
 EW-59981-21	1	0.063	230

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59981-22	1	0.053	270
 EW-59981-23	1	0.044	325

Cole-Parmer® 8 inch Dia. Sieves with NIST-Traceable Calibration Report



COLE-PARMER

Sieve openings are measured with a computerized video imaging system

Sieve is serialized, comes with a test sieve certificate, and meets ASTM E11 specifications. Includes an NIST-traceable calibration report supplied by the manufacturer.

item	Height (InchShort)	Material	U.S. Standard Screen Size (No.)
EW-59980-50	2	Brass frame, brass wire	8

item	Height (InchShort)	Material	U.S. Standard Screen Size (No.)
 EW-59980-51	2	Brass frame, brass wire	10
 EW-59980-52	2	Brass frame, brass wire	12
 EW-59980-53	2	Brass frame, brass wire	14

item	Height (InchShort)	Material	U.S. Standard Screen Size (No.)
 EW-59980-54	2	Brass frame, brass wire	16
 EW-59980-55	2	Brass frame, brass wire	18
 EW-59980-56	2	Brass frame, brass wire	20

item	Height (InchShort)	Material	U.S. Standard Screen Size (No.)
 EW-59980-57	2	Brass frame, brass wire	25
 EW-59980-58	2	Brass frame, brass wire	30
 EW-59980-59	2	Brass frame, brass wire	35

item	Height (InchShort)	Material	U.S. Standard Screen Size (No.)
 EW-59980-60	2	Brass frame, brass wire	40
 EW-59980-61	2	Brass frame, brass wire	45
 EW-59980-62	2	Brass frame, brass wire	50

item	Height (InchShort)	Material	U.S. Standard Screen Size (No.)
 EW-59980-63	2	Brass frame, brass wire	60
 EW-59980-64	2	Brass frame, brass wire	70
 EW-59980-65	2	Brass frame, brass wire	80

item	Height (InchShort)	Material	U.S. Standard Screen Size (No.)
 EW-59980-66	2	Brass frame, brass wire	100
 EW-59980-67	2	Brass frame, brass wire	120
 EW-59980-68	2	Brass frame, brass wire	140

item	Height (InchShort)	Material	U.S. Standard Screen Size (No.)
 EW-59980-69	2	Brass frame, brass wire	170
 EW-59980-70	2	Brass frame, brass wire	200
 EW-59980-71	2	Brass frame, brass wire	230

item	Height (InchShort)	Material	U.S. Standard Screen Size (No.)
 EW-59980-72	2	Brass frame, brass wire	270
 EW-59980-73	2	Brass frame, brass wire	325
 EW-59981-54	1	Brass frame, brass wire	16

item	Height (InchShort)	Material	U.S. Standard Screen Size (No.)
 EW-59981-58	1	Brass frame, brass wire	30
 EW-59981-59	1	Brass frame, brass wire	35
 EW-59981-60	1	Brass frame, brass wire	40

item	Height (InchShort)	Material	U.S. Standard Screen Size (No.)
 EW-59981-61	1	Brass frame, brass wire	45
 EW-59981-63	1	Brass frame, brass wire	60
 EW-59981-64	1	Brass frame, brass wire	70

item	Height (InchShort)	Material	U.S. Standard Screen Size (No.)
 EW-59981-65	1	Brass frame, brass wire	80
 EW-59981-66	1	Brass frame, brass wire	100
 EW-59981-73	1	Brass frame, brass wire	325

item	Height (InchShort)	Material	U.S. Standard Screen Size (No.)
 EW-59984-50	2	Stainless Steel	8
 EW-59984-51	2	Stainless Steel	10
 EW-59984-52	2	Stainless Steel	12

item	Height (InchShort)	Material	U.S. Standard Screen Size (No.)
 EW-59984-53	2	Stainless Steel	14
 EW-59984-54	2	Stainless Steel	16
 EW-59984-55	2	Stainless Steel	18

item	Height (InchShort)	Material	U.S. Standard Screen Size (No.)
 EW-59984-56	2	Stainless Steel	20
 EW-59984-57	2	Stainless Steel	25
 EW-59984-58	2	Stainless Steel	30

item	Height (InchShort)	Material	U.S. Standard Screen Size (No.)
 EW-59984-59	2	Stainless Steel	35
 EW-59984-60	2	Stainless Steel	40
 EW-59984-61	2	Stainless Steel	45

item	Height (InchShort)	Material	U.S. Standard Screen Size (No.)
 EW-59984-62	2	Stainless Steel	50
 EW-59984-63	2	Stainless Steel	60
 EW-59984-64	2	Stainless Steel	70

item	Height (InchShort)	Material	U.S. Standard Screen Size (No.)
 EW-59984-65	2	Stainless Steel	80
 EW-59984-66	2	Stainless Steel	100
 EW-59984-67	2	Stainless Steel	120

item	Height (InchShort)	Material	U.S. Standard Screen Size (No.)
 EW-59984-68	2	Stainless Steel	140
 EW-59984-69	2	Stainless Steel	170
 EW-59984-70	2	Stainless Steel	200

item	Height (InchShort)	Material	U.S. Standard Screen Size (No.)
 EW-59984-71	2	Stainless Steel	230
 EW-59984-72	2	Stainless Steel	270
 EW-59984-73	2	Stainless Steel	325

item	Height (InchShort)	Material	U.S. Standard Screen Size (No.)
 EW-59984-74	2	Stainless Steel	400
 EW-59985-50	1	Stainless Steel	8
 EW-59985-51	1	Stainless Steel	10

item	Height (InchShort)	Material	U.S. Standard Screen Size (No.)
 EW-59985-52	1	Stainless Steel	12
 EW-59985-54	1	Stainless Steel	16
 EW-59985-55	1	Stainless Steel	18

item		Height (InchShort)	Material	U.S. Standard Screen Size (No.)
 EW-59985-56		1	Stainless Steel	20
 EW-59985-57		1	Stainless Steel	25
 EW-59985-58		1	Stainless Steel	30

item	Height (InchShort)	Material	U.S. Standard Screen Size (No.)
 EW-59985-59	1	Stainless Steel	35
 EW-59985-60	1	Stainless Steel	40
 EW-59985-61	1	Stainless Steel	45

item	Height (InchShort)	Material	U.S. Standard Screen Size (No.)
 EW-59985-62	1	Stainless Steel	50
 EW-59985-63	1	Stainless Steel	60
 EW-59985-64	1	Stainless Steel	70

item	Height (InchShort)	Material	U.S. Standard Screen Size (No.)
 EW-59985-65	1	Stainless Steel	80
 EW-59985-66	1	Stainless Steel	100
 EW-59985-67	1	Stainless Steel	120

item	Height (InchShort)	Material	U.S. Standard Screen Size (No.)
 EW-59985-68	1	Stainless Steel	140
 EW-59985-69	1	Stainless Steel	170
 EW-59985-70	1	Stainless Steel	200

item	Height (InchShort)	Material	U.S. Standard Screen Size (No.)
 EW-59985-71	1	Stainless Steel	230
 EW-59985-72	1	Stainless Steel	270
 EW-59985-73	1	Stainless Steel	325

item	Height (InchShort)	Material	U.S. Standard Screen Size (No.)
 EW-59985-74	1	Stainless Steel	400

Cole-Parmer® Sieve Covers and Receiving Pans for 3" Diameter Sieves



- Suitable for use with 3” diameter sieves
- Constructed of durable stainless steel
Use covers and receiving pans with 3” diameter sieves. Order sieves separately. Order a receiving pan for the bottom sieve to collect the finest particles.

2 PRODUCT OPTIONS

item	Description
 EW-59988-03	3" Diameter Stainless Steel Sieve Cover
 EW-59988-07	3" Diameter Stainless Steel Full Height Sieve Pan

Cole-Parmer® 12 inch Dia. Sieves with NIST-Traceable Calibration Report



- **COLE-PARMER**

Sieve is serialized, comes with a test sieve certificate, and meets ASTM E11 specifications. Includes an NIST-traceable calibration report supplied by the manufacturer.

item	Height (InchShort)	Material	U.S. Standard Screen Size (No.)
 EW-59990-52	3.25	Brass frame, brass wire	12

item	Height (InchShort)	Material	U.S. Standard Screen Size (No.)
 EW-59990-55	3.25	Brass frame, brass wire	18
 EW-59990-59	3.25	Brass frame, brass wire	35
 EW-59991-51	1.625	Brass frame, brass wire	10

item	Height (InchShort)	Material	U.S. Standard Screen Size (No.)
 EW-59994-51	3.25	Stainless Steel	10
 EW-59994-53	3.25	Stainless Steel	14
 EW-59994-54	3.25	Stainless Steel	16

item	Height (InchShort)	Material	U.S. Standard Screen Size (No.)
 EW-59994-55	3.25	Stainless Steel	18
 EW-59994-56	3.25	Stainless Steel	20
 EW-59994-57	3.25	Stainless Steel	25

item	Height (InchShort)	Material	U.S. Standard Screen Size (No.)
 EW-59994-58	3.25	Stainless Steel	30
 EW-59994-59	3.25	Stainless Steel	35
 EW-59994-60	3.25	Stainless Steel	40

item	Height (InchShort)	Material	U.S. Standard Screen Size (No.)
 EW-59994-61	3.25	Stainless Steel	45
 EW-59994-62	3.25	Stainless Steel	50
 EW-59994-63	3.25	Stainless Steel	60

item	Height (InchShort)	Material	U.S. Standard Screen Size (No.)
 EW-59994-64	3.25	Stainless Steel	70
 EW-59994-65	3.25	Stainless Steel	80
 EW-59994-66	3.25	Stainless Steel	100

item	Height (InchShort)	Material	U.S. Standard Screen Size (No.)
 EW-59994-67	3.25	Stainless Steel	120
 EW-59994-68	3.25	Stainless Steel	140
 EW-59994-70	3.25	Stainless Steel	200

item	Height (InchShort)	Material	U.S. Standard Screen Size (No.)
 EW-59994-71	3.25	Stainless Steel	230
 EW-59994-73	3.25	Stainless Steel	325
 EW-59994-74	3.25	Stainless Steel	400

item	Height (InchShort)	Material	U.S. Standard Screen Size (No.)
 EW-59995-58	1.625	Stainless Steel	30

Cole-Parmer® Testing Sieves, 12" dia, Brass Frame/Brass Wire



These brass testing sieves are serialized and come with a test sieve certificate of compliance. Equivalent to ASTM E11-09 compliance sieve.

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59990-00	3.25	2.36	8

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59990-01	3.25	2	10
 EW-59990-02	3.25	1.68	12
 EW-59990-03	3.25	1.41	14

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59990-04	3.25	1.19	16
 EW-59990-05	3.25	1	18
 EW-59990-06	3.25	0.841	20

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59990-07	3.25	0.707	25
 EW-59990-08	3.25	0.595	30
 EW-59990-09	3.25	0.5	35

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59990-10	3.25	0.425	40
 EW-59990-11	3.25	0.355	45
 EW-59990-12	3.25	0.297	50

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59990-13	3.25	0.25	60
 EW-59990-14	3.25	0.212	70
 EW-59990-16	3.25	0.149	100

item		Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59990-18		3.25	0.105	140
 EW-59990-19		3.25	0.088	170
 EW-59990-20		3.25	0.074	200

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59990-23	3.25	0.044	325
 EW-59991-00	1.625	2.36	8
 EW-59991-01	1.625	2	10

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59991-02	1.625	1.68	12
 EW-59991-04	1.625	1.19	16
 EW-59991-05	1.625	1	18

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59991-09	1.625	0.5	35
 EW-59991-12	1.625	0.297	50
 EW-59991-14	1.625	0.212	70

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59991-16	1.625	0.149	100
 EW-59991-23	1.625	0.044	325

Cole-Parmer Testing Sieve, Brass, 200 mm OD, 0.500 mm Openings, Full Height, No. 35 mesh



image only.

ADD TO

Cole-Parmer – Item # EW-59942-18

Sieves are serialized and come with a test sieve certificate

- 1 year warranty

Specifications & Description

• Material	Brass frame, stainless steel wire
• U.S. Standard Screen Size (No.)	35
• Tyler Standard Screen Size (mesh)	32
• Sieve Opening (mm)	0.5
• Outside Diameter (in)	8
• Height (Inch/Short)	2
• Description	Testing Sieve, brass, 200 mm OD, 0.500 mm openings, full height, No. 35 mesh
• Warranty	1 Year



MORE ABOUT THIS ITEM

These 200-mm diameter sieves are the international equivalent to 8" diameter sieves. Use for dry particle analysis. Sieves feature stainless steel (SS) wire mesh and a sturdy brass frame. All sieves meet ISO 3310-1 and BS410 specifications. Stack these sieves to analyze several sizes of particles simultaneously.

Cole-Parmer® Full Height Sieves, Stainless Steel, 3" dia



Made of stainless steel (SS)

- Individually serial numbered
 - 1" stacking height
 - Compliant with ASTM E11 requirements for compliance sieves
- These test sieves feature a stainless steel woven-wire cloth available with various openings and a full height stainless steel frame. Ideal for use with Cole-Parmer sieve shaker (59986-01) , sold separately.

Optional accessories include a stainless steel cover (59988-03) and full-height receiving pan (59988-07) , sold separately.

•

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
EW-59987-01	1	4	5

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59987-02	1	3.35	6
 EW-59987-03	1	2.83	7
 EW-59987-04	1	2.36	8
 EW-59987-05	1	2	10

item		Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59987-06		1	1.68	12
 EW-59987-07		1	1.41	14
 EW-59987-08		1	1.19	16
 EW-59987-09		1	1	18

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59987-11	1	0.841	20
 EW-59987-12	1	0.707	25
 EW-59987-13	1	0.595	30
 EW-59987-14	1	0.5	35

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59987-16	1	0.425	40
 EW-59987-17	1	0.355	45
 EW-59987-18	1	0.297	50
 EW-59987-19	1	0.25	60

item		Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59987-20		1	0.212	70
 EW-59987-21		1	0.177	80
 EW-59987-22		1	0.149	100
 EW-59987-23		1	0.125	120

item	Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59987-24	1	0.105	140
 EW-59987-25	1	0.088	170
 EW-59987-26	1	0.074	200
 EW-59987-27	1	0.063	230

item		Height (InchShort)	Sieve Opening (mm)	U.S. Standard Screen Size (No.)
 EW-59987-28		1	0.053	270
		1	0.044	325
		1	0.037	400
EW-59987-29				
EW-59987-30				



По вопросам продаж и поддержки обращайтесь:

Алматы (727)345-47-04
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижегород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Саранск (8342)22-96-24
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(727)345-47-04

Беларусь +(375)257-127-884

Узбекистан +998(71)205-18-59

Киргизия +996(312)96-26-47

эл.почта: cen@nt-rt.ru || сайт: <http://coleparmer.nt-rt.ru/>